

DATASHEET

The Kim Infant CPR Manikin is a foam-filled and rugged manikin suitable for BLS training of a newborn, featuring anatomical landmarks including nipples, xiphoid process and sub-sternal notch.

The uni-body construction makes the manikin virtually indestructible and cross-contamination is controlled with the sanitary system.

Sanitary System

The individual mouth/nosepieces are easy to clean and the disposable airway system with a non-re-breathing valve prevents cross-contamination and is changed once a day.

Kim™ Infant CPR Manikin



KEY BENEFITS

- CPR training of newborns
- Lightweight and economical
- Tough and durable
- Anatomically correct with palpable landmarks
- Visible chest rise during ventilation
- Sanitary mouth/lung protection bags

ORDERING SPECIFICATIONS

Item no.	
100-2140U	Kim Infant CPR Manikin Full-body manikin with romper, 3 mouth/nosepieces, 2 airway systems, carry bag
100-2068U	Mouth/Nosepieces Kim (set of 10)
100-2905U	Kim/Kevin Airway System (set of 24)
100-2252U	Carry Bag for Kim Infant Manikin

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com